



GCMTI RD-2:2024

**Determination of Psoralen and Isopsoralen
in Proprietary Chinese Medicines containing Psoraleae and Ginseng
by Liquid Chromatograph-Tandem Mass Spectrometer
(LC-MS/MS)**

GCMTI Method Publications



Determination of Psoralen and Isopsoralen
in Proprietary Chinese Medicines containing Psoraleae and Ginseng
by Liquid Chromatograph-Tandem Mass Spectrometer
(LC-MS/MS)¹

***Safety Precaution:** This procedure involves carcinogenic chemicals, corrosive chemicals and flammable solvents. Apply precautions when handling such chemicals, for example: use eye and hand protection and where necessary carry out the work in a fume cupboard to avoid inhalation of vapour.*

1. Introduction

- 1.1. Proprietary Chinese medicines (pCm) containing psoraleae (補骨脂) and ginseng (人參) for nourishing and Yang invigorating are commonly found in Hong Kong. Nevertheless, analysis of the chemical markers of psoraleae and ginseng is a great challenge since it is susceptible to interference from matrix and other chemical components.
- 1.2. This method describes the procedures for qualitative and quantitative determination of psoralen and isopsoralen in pCm containing psoraleae and ginseng by liquid chromatograph-tandem mass spectrometer (LC-MS/MS).

2. Reagents

Note: All reagents used should be of analytical reagent grade or equivalent unless otherwise specified.

- 2.1. Methanol, LC-MS grade.
- 2.2. Milli-Q water.
- 2.3. Psoralen, CAS no.: 66-97-7.
- 2.4. Isopsoralen, CAS no.: 523-50-2.
- 2.5. Extraction solvent

Methanol : water (7:3 v/v).

¹ This method intends to provide a reliable analytical method that can be used as a quality control method for determining the targeted chemical marker(s) in the corresponding pCm product(s). It is the user's responsibility to assess the suitability of testing their pCm products when adopting this method.

2.6. Preparation of standard solutions

2.6.1. Individual stock standard solutions (ca. 1000 µg/mL)

Weigh accurately about 10 mg of psoralen (Clause 2.3) and isopsoralen (Clause 2.4) into separate 10-mL volumetric flasks, dissolve and make up to the mark with methanol (Clause 2.1), respectively.

2.6.2. Mixed intermediate standard solution I (ca. 30 µg/mL)

Transfer 0.3 mL of each individual stock standard solutions into a 10-mL volumetric flask and make up to the mark with methanol (Clause 2.1).

2.6.3. Mixed intermediate standard solution II (ca. 120 ng/mL)

Transfer 0.1 mL of mixed intermediate standard solution I into a 25-mL volumetric flask and make up to the mark with extraction solvent (Clause 2.5).

2.6.4. Calibration standard solutions, CS1 – CS5

A series of calibration standard solutions are prepared by transferring an appropriate amount of mixed intermediate standard solution II into 10-mL volumetric flasks and make up to the mark with extraction solvent (Clause 2.5). Suggested volumes of standard solution used for the preparation are listed in the table below.

Calibration standard	Volume of mixed intermediate standard solution II (mL)	Final Volume (mL)	Conc. of psoralen and isopsoralen (ng/mL)
CS1	0.25	10	3
CS2	0.50	10	6
CS3	0.75	10	9
CS4	1.00	10	12
CS5	1.25	10	15

2.6.5. Individual stock initial calibration verification (ICV) standard solutions (ca. 1000 µg/mL)

Prepare individual stock ICV standard solutions, from source different from that of the calibration standard. Weigh accurately about 10 mg of psoralen (Clause 2.3) and isopsoralen (Clause 2.4) into separate 10-mL volumetric flasks, dissolve and make up to the mark with methanol (Clause 2.1), respectively.

2.6.6. Mixed intermediate ICV standard solution I (ca. 30 µg/mL)

Transfer 0.3 mL of each individual stock ICV standard solutions into a 25-mL volumetric flask and make up to the mark with methanol (Clause 2.1).

2.6.7. Mixed intermediate ICV standard solution II (ca. 120 ng/mL)

Transfer 0.1 mL of mixed intermediate ICV standard solution I into a 25-mL volumetric flask and make up to the mark with extraction solvent (Clause 2.5).

2.6.8. ICV working standard solution (ca. 9 ng/mL)

Transfer 0.75 mL of mixed intermediate ICV standard solution III into a 10-mL volumetric flask and make up to the mark with extraction solvent (Clause 2.5).

2.6.9. Spike standard solutions (ca. 1000 µg/mL)

Refer to individual stock standard solutions (Clause 2.6.1).

3. Apparatus

All glassware shall be rinsed with acetone and washed with detergent solution as soon as practicable after use. After detergent washing, glassware shall be rinsed immediately, firstly with water and then with acetone twice.

3.1. Grinder or blender.

3.2. Analytical balance, capable of weighing to 0.01 mg.

3.3. Volumetric flasks, 10-mL and 25-mL.

3.4. Auto pipettes, 100-µL, 300-µL and 1000-µL.

3.5. Centrifuge with rotation speed of at least 4000 rpm.

3.6. Centrifuge tubes, 15-mL.

3.7. Vortex mixer.

3.8. Ultrasonic bath.

3.9. PTFE membrane filters, 0.2 µm.

3.10. LC glass vials.

3.11. LC column: InertSustain C18, 5 µm, 2.1 mm × 250 mm, GL Sciences or equivalent.

3.12. Liquid chromatograph-tandem mass spectrometer (LC-MS/MS) system.

4. Procedures

4.1. Sample extraction

- 4.1.1. Grind and homogenise solid samples using grinder or blender.
- 4.1.2. Weigh accurately about 0.25 g of sample into a 15-mL centrifuge tube.
- 4.1.3. Add 10 mL of extraction solvent (Clause 2.5) into the centrifuge tube. Vortex the sample mixture for 1 minute.
- 4.1.4. Sonicate the sample mixture in an ultrasonic bath for 20 minutes at room temperature.
- 4.1.5. Centrifuge the sample solution at 4000 rpm for 10 minutes. Carefully transfer the supernatant solution to a 25-mL volumetric flask.
- 4.1.6. Repeat Clauses 4.1.3 to 4.1.5 twice with 5 mL of extraction solvent (Clause 2.5). Collect all supernatant in the same 25-mL volumetric flask and make up to the mark with extraction solvent (Clause 2.5). Dilute the sample solution by 50-fold with extraction solvent (Clause 2.5).
- 4.1.7. Filter the diluted sample solution with 0.2 µm PTFE membrane filter into a LC glass vial. The solution is ready for LC-MS/MS analysis.

Remark:

Dilute the sample solution with extraction solvent (Clause 2.5) if the concentration of analyte(s) is not within the calibration range.

4.2. LC-MS/MS analysis

- 4.2.1. Operate the LC-MS/MS system in accordance with the instrument manual. Carry out analysis with the conditions as suggested below. It may be necessary to modify the operation conditions for optimal signal output. Record the actual experimental conditions in the worksheet.

4.2.2. Suggested LC conditions:

LC system	:	Thermo Scientific UltiMate 3000 UHPLC or equivalent performance
Column	:	GL Sciences InertSustain C18, 5 µm, 2.1 mm × 250 mm or equivalent
Column temperature	:	40 °C
Flow rate	:	0.3 mL/min
Injection volume	:	5 µL
Mobile phase	:	A: Water B: Methanol

Gradient	:	Time (min)	A%	B%
		0.0	60	40
		2.0	60	40
		16.0	30	70
		18.0	30	70
		19.0	15	85
		22.0	15	85
		22.1	5	95
		25.0	5	95
		25.1	60	40
		28.0	60	40

4.2.3. Suggested MS/MS conditions:

MS/MS system	:	SCIEX 6500+ system
Ionization mode	:	Electrospray ionization (ESI)
Polarity	:	Positive mode
Ionspray voltage	:	4500V
Source temperature	:	350°C
Ion source gas 1 (GS1)	:	40
Ion source gas 2 (GS2)	:	40
Curtain gas (CUR)	:	25
Collision gas (CAD)	:	Medium
Scan type	:	MRM

4.2.4. Suggested MRM acquisition parameters:

Analyte	MRM transition	Dwell time (msec)	DP	EP	CE	CXP
Psoralen	186.9 → 131.0*	50	91	10	33	26
	186.9 → 115.0 [^]	50	91	10	31	16
Isopsoralen	186.9 → 131.0*	50	91	10	31	14
	186.9 → 115.0 [^]	50	91	10	31	16

*Remark: The quantification MRMs and the qualification MRMs are marked with * and [^] respectively.*

4.2.5. Calibrate the LC-MS/MS system using at least 5 calibration standards (Clause 2.6.4).

4.2.6. Perform LC-MS/MS analysis for method blank(s), sample(s), sample duplicate(s), spike sample(s) and relevant check standard solution(s) according to the quality control plan as established in the laboratory.

5. Calculation / result interpretation

5.1. Identification requirement

5.1.1. Identify the target analyte in the sample by comparison of the retention

time of the detected peak (RT_{sample}) with that of the average retention time (RT) of the calibration standards. The RT_{sample} shall not differ from that of the average RT of calibration standards by more than 5% for positive identification.

- 5.1.2. The relative abundance of MRMs shall meet the tolerance for positive identification of the analyte(s) (with reference to that of the average relative abundance of the calibration standard):

Relative intensity to the base peak	% Allowable deviation
>50%	±20%
>20% to 50%	±25%
>10% to 20%	±30%
≤10%	±50%

- 5.2. Establish the calibration curve by plotting the peak area against the concentration of analyte in the calibration standards in linear calibration mode.
- 5.3. Calculate the concentration of analyte in the sample, in $\mu\text{g/g}$, using the following equation:

$$\text{Concentration of analyte } (\mu\text{g/g}) = \frac{C \times V \times D}{1000 \times W}$$

where C = Conc. of analyte obtained from calibration curve (in ng/mL)
 V = Final volume (mL)
 D = Dilution factor
 W = Sample weight (g)

- 5.4. If matrix effect is suspected when significant bias is detected in spike recovery, it may be minimized by (1) further dilution of the sample solution or (2) quantification using standard addition approach.

6. Reference

- 6.1. Chinese Pharmacopoeia Commission. Pharmacopoeia of the People's Republic of China Volume 1, 2020 ed. China Medical Science Press.
- 6.2. "Quantifying Uncertainty in Analytical Measurement", Eurachem/ CITAC Guide CG4, 3rd Edition, 2012.
- 6.3. V. J. Barwick and S. L. R. Ellison, "VAM Project 3.2.1 Development and Harmonisation of Measurement Uncertainty Principles Part (d): Protocol for Uncertainty Evaluation from Validation data", LGC/VAM/1998/088 Version 5.1, January 2000.